



## TECHNICAL DATA SHEET

www.micromod.de  
Telephone: + 49 381 / 7999 7000  
Fax: + 49 381 / 7999 7070  
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH  
Schillingallee 68  
D-18057 Rostock  
Germany

---

<b>Product code:</b>	01-05-303
<b>Product name:</b>	micromer®
<b>Surface:</b>	NR3+
<b>Size:</b>	3 µm
<b>Solid content:</b>	50 mg/ml
<b>Composition:</b>	polystyrene matrix particles
<b>Standard deviation:</b>	< 5 % (C.V.)
<b>Shape:</b>	spherical
<b>Density:</b>	1.03 g/ccm
<b>Stable in:</b>	aqueous buffers, methanol, ethanol, DMSO
<b>Not stable in:</b>	halogenated hydrocarbons, toluene, strong acidic solutions, e.g. 10% HCl
<b>Product form:</b>	Suspension in water
<b>Particles per ml:</b>	3.4 E9
<b>Particles per mg:</b>	6.9 E7
<b>Colour:</b>	white
<b>Additional remarks:</b>	Storage at room temperature.