



TECHNICAL DATA SHEET

www.micromod.de
Telephone: + 49 381 / 7999 7000
Fax: + 49 381 / 7999 7070
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
Schillingallee 68
D-18057 Rostock
Germany

| | |
|----------------------------|---|
| Product code: | 01-25-303 |
| Product name: | micromer® |
| Surface: | collagen |
| Size: | 3 µm |
| Solid content: | 25 mg/ml |
| Composition: | polystyrene matrix particles |
| Standard deviation: | < 5 % (C.V.) |
| Shape: | spherical |
| Density: | 1.03 g/ccm |
| Stable in: | aqueous buffers pH > 4 |
| Not stable in: | organic solvents, acidic solutions pH < 4 |
| Product form: | Suspension in water |
| Particles per ml: | 1.7 E9 |
| Particles per mg: | 6.9 E7 |
| Colour: | white |
| Additional remarks: | Storage at 2 - 8 °C for 6 months. |