



TECHNICAL DATA SHEET

www.micromod.de
Telephone: + 49 381 / 7999 7000
Fax: + 49 381 / 7999 7070
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
Schillingallee 68
D-18057 Rostock
Germany

Product code:	104-19-501
Product name:	synomag®-D
Surface:	streptavidin
Size:	50 nm
Solid content:	5 mg/ml
Composition:	dextran iron oxide composite particles
Polydispersity index:	< 0.200
Shape:	cluster-typed
Density:	2.5 mg/ml
Iron concentration:	3 mg/ml
Magnetization:	100 Am ² /kg iron (H = 80 kA/m)
Saturation magnetization:	> 110 Am ² /kg iron (H = 800 kA/m)
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	Suspension in PBS (0.02 % sodium azide)
Particles per ml:	3.1 E13
Particles per mg:	6.1 E12
Biotin binding capacity:	> 200 pmol/mg iron
Colour:	dark brown
Additional remarks:	Storage at 2 - 8 °C for 3 months.