



TECHNICAL DATA SHEET

www.micromod.de
Telephone: + 49 381 / 7999 7000
Fax: + 49 381 / 7999 7070
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
Schillingallee 68
D-18057 Rostock
Germany

Product code:	11-00-502
Product name:	PLA particles
Surface:	plain
Size:	500 nm
Solid content:	10 mg/ml
Composition:	poly(lactic acid) particles
Polydispersity index:	< 0.200
Shape:	spherical
Density:	1.0 g/ccm
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	Suspension in water
Particles per ml:	1.5 E11
Particles per mg:	1.5 E10
Colour:	white
Additional remarks:	Storage at 2 - 8 °C for 6 months.