



TECHNICAL DATA SHEET

www.micromod.de
Telephone: + 49 381 / 7999 7000
Fax: + 49 381 / 7999 7070
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
Schillingallee 68
D-18057 Rostock
Germany

Product code:	126-01-701
Product name:	synomag®-CLD-far redF
Surface:	NH ₂
Size:	70 nm
Solid content:	5 mg/ml
Composition:	cross-linked dextran iron oxide composite particles
Polydispersity index:	< 0.200
Shape:	cluster-typed
Density:	2.5 g/ccm
Iron concentration:	3 mg/ml
Magnetization:	100 Am ² /kg iron (H = 80 kA/m)
Saturation magnetization:	> 110 Am ² /kg iron (H > 800 kA/m)
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	Suspension in PBS buffer (pH = 7.4)
Particles per ml:	1.1 E13
Particles per mg:	2.2 E12
Optical properties:	far-red fluorescent
Colour:	dark brown
Excitation:	732 nm
Emission:	758 nm
Additional remarks:	Storage at 2 - 8 °C for 6 months.
TECHNOTES:	The particle matrix consists of cross-linked dextran.