



TECHNICAL DATA SHEET

www.micromod.de
 Telephone: + 49 381 / 7999 7000
 Fax: + 49 381 / 7999 7070
 E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
 Schillingallee 68
 D-18057 Rostock
 Germany

Product code:	194-19-102
Product name:	BNF-Dextran-far redF
Surface:	streptavidin
Size:	100 nm
Solid content:	5 mg/ml
Composition:	cross-linked dextran iron oxide composite particles
Polydispersity index:	< 0.250
Shape:	cluster-typed
Density:	3.2 g/ccm
Iron concentration:	3 mg/ml
Magnetization:	49 Am ² /kg iron (H = 80 kA/m)
Saturation magnetization:	>76 Am ² /kg iron (H > 800 kA/m)
Coercive field H_c:	0.449 kA/m
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	Suspension in PBS (0.02 % sodium azide)
Particles per ml:	3.0 E12
Particles per mg:	6.0 E11
Biotin binding capacity:	> 100 pmol/mg iron
Optical properties:	far-red fluorescent
Colour:	dark brown
Excitation:	732 nm
Emission:	758 nm
Additional remarks:	Storage at 2 - 8 °C for 3 months.
TECHNOTES:	The particle matrix consists of crosslinked dextran.