



TECHNICAL DATA SHEET

www.micromod.de
Telephone: + 49 381 / 7999 7000
Fax: + 49 381 / 7999 7070
E-mail address: info@micromod.de

micromod Partikeltechnologie GmbH
Schillingallee 68
D-18057 Rostock
Germany

Product code:	51-02-105
Product name:	PLA-greenF
Surface:	COOH
Size:	100 µm
Solid content:	10 mg/ml
Composition:	green fluorescent poly(lactic acid) particles
Standard deviation:)*
Shape:	spherical
Density:	1.0 g/ccm
Stable in:	aqueous buffers pH > 4
Not stable in:	organic solvents, acidic solutions pH < 4
Product form:	Suspension in water
Particles per ml:	1.9 E4
Particles per mg:	1.9 E3
Optical properties:	green fluorescent
Colour:	pink
Excitation:	502 nm
Emission:	527 nm
Additional remarks:	Storage at 2 - 8 °C for 6 months.
TECHNOTES:)* Size cut: 70 µm - 150 µm